homecare













JOYCEeasy












Nasal Mask

Full Face Mask

Instructions for Use

WEINMANN
medical technology

Deutsch	3
English	14
Français	24
Nederlands	35
Italiano	46
Español	57
Português	68
Русский	79
Polski	91
Hrvatski	102
Türkçe	113
Ελληνικά	123
中文	135

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Table of Contents

1 Introduction	15	9 Technical data	21
1.1 Intended use	15	9.1 Materials	22
1.2 Contraindications	15	9.2 Characteristic pressure/ flow curve	23
1.3 Side effects	15	10 Warranty	23
2 Safety	15	11 Declaration of conformity	23
2.1 Safety information	15		
2.2 General information	16		
3 Product description	16		
3.1 Overview	16		
3.2 Compatible devices	16		
3.3 Exhalation system	17		
3.4 Anti-asphyxia valve	17		
3.5 Pressure measurement port	18		
4 Operation	18		
4.1 Applying the mask	18		
4.2 Adjusting the mask	18		
4.3 Removing the mask	18		
4.4 Dismantling the mask	18		
4.5 Assembling the mask	19		
5 Hygiene preparation	19		
5.1 Cleaning the mask	19		
5.2 Tips	20		
5.3 Change of patient	20		
6 Disposal	20		
7 Faults	21		
8 Scope of supply	21		

1 Introduction

1.1 Intended use

JOYCEeasy nasal masks and JOYCEeasy Full Face masks are used for treating sleep apnea and for non-invasive and non-life-sustaining ventilation of patients with respiratory insufficiency. It serves as a connecting element between the patient and the therapy device.

1.2 Contraindications

In the case that any of the following symptoms are present the mask may not be used, or only with particular caution:

erosions and ulcerations, skin allergies, flush of the facial skin, pressure marks on the face, claustrophobia, anxiety, deformities of the face or nasopharynx, intake of medication which may induce vomiting, the necessity for immediate intubation. Please also observe the contraindications in the instructions for use of your therapy device.

1.3 Side effects

Nasal congestion, dry nose, morning xerostomia (dry mouth), feeling of pressure in the sinuses, irritated conjunctiva, skin rashes, pressure marks in the face, irritating noises when breathing.

2 Safety

Please read these instructions carefully. They are a constituent part of the mask and must be available at all times. In accordance with Directive 93/42/EEC, observe the points below.

2.1 Safety information

Risk of injury from the supply of oxygen!

Oxygen can deposit on clothing, bed linen and hair. In conjunction with smoking, naked flames or electrical equipment, it can cause fires and explosions.

- ⇒ Do not smoke.
- ⇒ Avoid naked flames.
- ⇒ Do not use electrically conductive tubes (exception: breathing tubes which can be heated electrically).

Risk of injury from oxygen therapy!

Oxygen therapies can lead to side effects in case of incorrect use or dosage.

- ⇒ Only administer oxygen as prescribed by a physician.
- ⇒ Follow the instructions for use of the oxygen system and therapy device.

Risk of injury from re-inhaling CO₂!

If the mask is used incorrectly, CO₂ may be re-inhaled.

- ⇒ Do not close the openings of the mask.
- ⇒ Only apply the mask for a prolonged period if the therapy device is running.
- ⇒ Only use the mask within the quoted therapy pressure range.
- ⇒ Do not use the mask on patients who cannot take it off themselves.

⇒ Monitor patients with obstructive and restrictive lung diseases individually when using JOYCEeasy masks.

Risk of injury from slipping of the mask!

If the mask slips or falls off, the therapy is ineffective.

- ⇒ Monitor patients with restricted spontaneous respiration.
- ⇒ Activate low pressure/leakage alarms on the therapy device.
- ⇒ Monitor patients with obstructive and restrictive lung diseases individually.

Patient endangerment through application of mask during anesthesia!

When the nasal mask is used during anesthesia, the anesthetic gas supplied can escape through the mouth and endanger the patient.

- ⇒ Do not use nasal mask during anesthesia.

2.2 General information

- Use of third-party products may lead to functional failures and restricted fitness for purpose. Biocompatibility may also be compromised. Please note that in these cases, any claim under warranty and liability will be void if neither the accessories nor original replacement parts recommended in the instructions are used.
- To prevent infection, bacterial contamination or functional impairments, follow the instructions in the section "Hygiene preparation" (see "5 Hygiene preparation", page 19).

3 Product description

3.1 Overview

Illustrations of the individual components can be found on the enclosed folding sheet (**figure A1**).

Key

- 1** Headgear
- 2** Forehead cushion
- 3** Forehead support
- 4** Lock
- 5** Retaining ring
- 6** Elbow
- 7** Plug for pressure measurement port
- 8** Pressure measurement port
- 9** Rotating sleeve
- 10** Mask body
- 11** Headgear clip
- 12** Mask cushion
- 13** Anti-asphyxia valve (on full-face masks only, see **figure F**)

3.2 Compatible devices

The mask can be used with all therapy devices which are not for life-sustaining treatment of patients with sleep apnea or ventilatory insufficiency.

With many combinations of device, the actual pressure in the mask will not correspond to the prescribed therapy pressure, even if the therapy device is displaying the correct pressure. Have the combination of devices adjusted by a physician or specialist dealer so that the actual pressure in the mask corresponds to therapy pressure.

Quick-release cord (optional)

You can use a quick-release cord (accessory) to open the headgear in emergency situations (e.g. vomiting).

The quick-release cord consists of a special clip with a cord which you can fit to the mask instead of one of the bottom headgear clips.

In an emergency, simply pull the quick-release cord. The special clip will detach itself from the mask. The headgear is now open and the mask can be removed (**figure G**).

3.3 Exhalation system

The mask has an integrated exhalation system. The retaining ring and mask body are shaped in such a way that a gap forms between these parts. The exhaled air can escape through this gap (**figure A2**).

3.4 Anti-asphyxia valve

WARNING

Risk of asphyxia if anti-asphyxia valve not working properly!

Residues may block the mask and jeopardize the success of therapy.

- ⇒ Check before each use that the anti-asphyxia valve is working properly.
- ⇒ Clean contact surfaces between the elbow and the rotating sleeve thoroughly.

Safety function of the anti-asphyxia valve

The anti-asphyxia valve has an opening to the atmosphere through which the patient can inhale ambient air if the therapy device fails. This significantly reduces the risk of breathing CO₂ back in and thus also the risk of asphyxia. Inside the anti-asphyxia valve is a membrane which can adopt one of two positions:

- as long as the therapy device is in operation, the air flow forces the membrane against the opening to the atmosphere and seals it completely. The patient breathes through the therapy device.
- if the therapy device fails or if the device is not connected, the membrane is in the rest position, in other words the opening to the atmosphere is free. The patient breathes ambient air through this opening.

3.5 Pressure measurement port

The mask has a connection **8** which can be used for measuring pressure or supplying oxygen. If you are not using the pressure measurement port, close it with the plug **7** to reach therapy pressure.

4 Operation

4.1 Applying the mask

1. Pull the headgear **1** over the head and place the mask on the face (**figure B1**).
2. Insert headgear clips **11** (**figure B2**).

4.2 Adjusting the mask

1. Undo lock **4** (**figure C1**).
2. Pull the forehead support **3** in or out of the lock **4** bit by bit as necessary (**figure C1**).

Remember the most convenient position. To do this, use the scale between the forehead support and the lock (**figure C1**).

3. Applying the mask (see "4.1 Applying the mask", page 18).
4. Adjust headgear **1** with the aid of the hook-and-loop fasteners so that the mask seals on the face but is not too tight (**figure C2**).

5. To adjust the hook-and-loop fastener on top of the head undo the hook-and-loop fastener, pull on the ends of the headgear with both hands and close the fastener again (**figure C3**).

4.3 Removing the mask

1. Undo **11** the headgear clips.
2. Slip the headgear **1** over the head and remove the mask.

4.4 Dismantling the mask

1. Undo the headgear clips **11** (**figure D1**).
2. Release hook-and-loop fasteners from headgear clips **11** and forehead support **3** and remove the headgear **1** (**figure D1 + D2**).
3. Remove the plug for the pressure measurement port **7** (**figure D3**).
4. Remove forehead cushion **2** from forehead support **3** (**figure D4**).
5. Take mask cushion **12** off mask body **10** (**figure D5**).
6. Turn retaining ring **5** counterclockwise and remove it (**figure D6**).
7. Remove elbow **6** from retaining ring **5** (**figure D7**).
8. Take rotating sleeve **9** off elbow **6** (**figure D7**).

- For full-face masks only: carefully take the anti-asphyxia valve off the elbow. In the process, please note: do not damage the membrane on the inside of the anti-asphyxia valve.

4.5 Assembling the mask

- Push rotating sleeve **9** onto elbow **6** (**figure E1**).
- Insert elbow **6** in the retaining ring **5** (**figure E1**).
- Insert retaining ring in the mask body with position **5** at 11 o'clock and turn clockwise to a 12 o'clock position until it locks into place (**figure E2**).
- Fit mask cushion **12** onto mask body **10** (**figure E3**).
- Pull forehead cushion **2** over the forehead support **3** (**figure E4**).
- Fit the plug for the pressure measurement port **7** (**figure E5**).
- Attach the headgear **1** with the headgear clips using the lower hook-and-loop fasteners **11** (**figure E6**).
- Attach the headgear **1** to the mask (**figure E7**).
- For full-face masks only: carefully push anti-asphyxia valve onto the elbow until it audibly engages on both sides of the elbow (**figure F**).
In the process, please note: do not damage or trap the membrane on the inside of the anti-asphyxia valve.

You can tell whether the headgear is correctly attached by the hook-and-loop tabs. The hook-and-loop tabs on the forehead have the same color as the inside of the head bands.

5 Hygiene preparation

WARNING

Risk of injury from inadequate cleaning!

Residues may congest the mask, impair the integrated exhalation system and jeopardize therapy success.

⇒ For patients with a compromised immune system or particular background of illness, disinfect mask components daily following consultation with the physician.

⇒ Thoroughly clean the contact surfaces between the elbow and the rotating sleeve.

- This product may contain disposable items. Disposable items are intended for single use.** Therefore, use these items only once and do **not** reprocess them. Reprocessing disposable items may impair functionality and safety of the product and lead to unforeseeable reactions due to ageing, embrittlement, wear, thermal stress, the effects of chemical processes etc.

5.1 Cleaning the mask

- Dismantle the mask (see "4.4 Dismantling the mask", page 18).

2. Clean mask in accordance with the table below:

Action	Daily	Weekly
Wash mask components with warm water and a mild detergent.	X	
When washing mask parts (exception: mask cushion, anti-asphyxia valve and forehead cushion), clean thoroughly with a cloth or a soft brush.		X
Wash headgear by hand.		X

3. Rinse all parts with clear water.

NOTICE

Material damage due to abrasion!

Rubbing the mask cushion and forehead cushion can damage the coating.
 ⇒ Do not rub the mask cushion and forehead cushion dry.

NOTICE

Material damage due to careless use!

Careless handling can damage the anti-asphyxia valve.
 ⇒ Always handle the anti-asphyxia valve with care.

4. Allow all parts to air-dry.
5. Perform a visual inspection.
6. Replace damaged parts if necessary.
7. Assembling the mask (see "4.5 Assembling the mask", page 19).

i Discolorations of mask components do not impair the functionality of the mask.

5.2 Tips

- You can wash all parts (exception: anti-asphyxia valve and headgear) in a dishwasher at 65 °C.
- Do not dry parts in direct sunlight.
- Wash headgear before first use, as the color may run.
- Do not iron the headgear in order to keep the hook-and-loop fasteners intact.
- Do not dry the headgear in a tumble-dryer.

5.3 Change of patient

You must subject the mask to a hygiene preparation in the event of a change of patient. For information on hygiene preparation in the event of a change of patient, please refer to our internet brochure at www.weinmann.de. If you do not have Internet access, please contact Weinmann.

6 Disposal

You can dispose of all parts in domestic waste.

7 Faults

Fault	Cause	Remedy
Facial pressure pain.	Mask fits too tightly.	Loosen headgear. Adjust position of forehead support.
Draft in the eyes.	Mask fits too loosely.	Tighten headgear. Adjust position of forehead support.
	Mask does not fit.	Use different mask size. Contact your specialist retailer.
Therapy pressure is not reached.	Mask is not adjusted correctly.	Re-adjust mask (see 4.2, p. 18).
	Mask cushion is damaged.	Replace mask cushion.
	Tube system is leaking.	Check connection and fit of tubes.
	Air is escaping from pressure measurement connection.	Close pressure measurement port with plug (included in scope of supply).
	Anti-asphyxia valve incorrectly fitted.	Fit anti-asphyxia valve correctly (see 4.5, p. 19).
	Anti-asphyxia valve defective	Replace anti-asphyxia valve

8 Scope of supply

Size	Variant	Nasal cannula	Full-face mask
S	vented	WM 25510	WM 25910
M	vented	WM 25520	WM 25920
L	vented	WM 25530	WM 25930

You can order replacement parts separately if necessary. A current list of replacement parts can be ordered on the Internet at www.weinmann.de or via your specialist dealer.

9 Technical data

	Nasal cannula	Full-face mask
Product class according to directive 93/42/EEC	IIa	
Dimensions (W x H x D)	approx. 80 mm x 130 mm x 90 mm ¹	approx. 120 mm x 150 mm x 110 mm ¹
Weight	approx. 72 g ¹	approx. 103 g ¹
Therapy pressure range	4 hPa - 30 hPa	
Hose connection: Cone according to EN ISO 5356-1	Ø 22 mm (male)	
Pressure measurement port	Ø 4 mm	
Wide headgear clips	max. 20 mm	

	Nasal cannula	Full-face mask
Temperature range: Operation Storage	+5 °C to + 40 °C -20 °C to +70 °C	
Flow resistance at 50 l/min at 100 l/min	0.16 hPa 0.64 hPa	0.13 hPa 0.53 hPa
Anti-asphyxia valve flow resistance • Inspiration at 50 l/min: • Exhalation at 50 l/min:	-	0.5 hPa 0.5 hPa
Switching pressure for anti-asphyxia valve • Open: • Close:	-	≤ 1 hPa ≤ 2 hPa
Quoted two-figure noise emission value according to ISO 4871: -Sound pressure level -Sound power level, uncertainty factor:	18 dB(A) 23 dB(A) 3 dB(A)	
Service life	Up to 12 months ²	
Standards applied	EN ISO 17510-2:2009	

¹ depending on size and variant.

² The materials of the mask deteriorate if exposed to e.g. aggressive detergents. In individual cases, therefore, it may be necessary to replace mask parts sooner.

CE 0197 The right to make design modifications is reserved.


9.1 Materials

Mask part	Material
Mask cushion, forehead cushion	Silicone, silicone coating
Plug for pressure measurement port	Silicone
Retaining ring (hard component)	PP (polypropylene)
Retaining ring (soft component)	TPE (thermoplastic elastomer)
Headgear clip, lock	PA (polyamide), POM (polyoxymethylene)
Forehead support, rotating sleeve, mask body, elbow	PA (polyamide)
Headgear	Elastane, polyester, PU (polyurethane), UBL loop, cotton, PA (polyamide)
Anti-asphyxia valve	TPE, PP (polypropylene)
Quick-release cord	POM, polyester, PA (polyamide)

All parts of the mask do not contain latex, PVC (polyvinylchloride) and DEHP (diethylhexyl phthalate).

9.2 Characteristic pressure/flow curve

The characteristic pressure/flow curve shows the outlet flow in dependence with therapy pressure.

10 Warranty

Weinmann gives the customer a limited manufacturer warranty on new genuine Weinmann products and any replacement part fitted by Weinmann in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. The warranty conditions can be downloaded from www.weinmann.de on the Internet. We can also send you the warranty conditions on request. In the event of a claim under warranty, please contact your specialist dealer.

Product	Period of guarantee
Weinmann devices including accessories (except masks) for sleep diagnosis, sleep therapy, home ventilation, oxygen medicine and emergency medicine	2 years
Masks including accessories, rechargeable batteries, batteries*, sensors, tube systems	6 months
Disposable products	None

*unless there is information to the contrary in the technical documentation.

11 Declaration of conformity

Weinmann Geräte für Medizin GmbH + Co. KG hereby declares that the product complies with the relevant regulations of Directive 93/42/EEC governing medical devices. The complete text of the declaration of conformity can be found at: www.weinmann.de

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Geräte für Medizin GmbH+Co.KG